


BAYMAX : *Personal Medical Assistant*



SYSTEM ARCHITECTURE





BAYMAX



Home



Symptoms



Lab Report




X-Ray Report


Symptom-based Medical Prediction

Enter Symptoms (comma-separated):

Predict

 BAYMAX

 Home

 Symptoms

 **Lab Report**

 X-Ray Report

OCR Lab Report Analysis

Upload an image file to extract text using OCR.

Choose an image file...



Drag and drop file here

Limit 200MB per file • PNG, JPG, JPEG, BMP

Browse files



Screenshot 2024-09-26 172846.png 68.7KB



Yash M. Patel

Age : 21 Years

Sex : Male

PID : 555



Sample Collected At:

125, Shivam Bungalow, S G Road,
Mumbai

Ref. By: Dr. Hiren Shah



Registered on: 02:31 PM 02 Dec, 2X

Collected on: 03:11 PM 02 Dec, 2X

Reported on: 04:35 PM 02 Dec, 2X

PERIPHERAL BLOOD SMEAR EXAMINATION


Investigation


Result

Reference Value

Unit

 BAYMAX

 Home

 Symptoms

 Lab Report

 **X-Ray Report**

Pneumonia classification

Please upload a chest X-ray image



Drag and drop file here

Limit 200MB per file • JPEG, JPG, PNG

[Browse files](#)



NORMAL2-IM-1427-0001.jpeg 247.2KB



SYMPTOM-BASED PREDICTION

- The symptom-based prediction feature allows users to input their symptoms to receive potential diagnoses and tailored recommendations.
- It leverages machine learning to analyze input data and provide accurate results.

Enter Symptoms (comma-separated):

itchin,skinrasehs

Predict

Normalized Symptoms:

itching, skin_rash

Predicted Disease: Fungal infection

HOW IT WORKS :

Step 1: User inputs symptoms in a comma-separated format (e.g., fever, cough).

Step 2: The system normalizes symptoms through:

- Tokenization: Breaks symptoms into individual words.

- Stop Word Removal: Filters out common words (e.g., "and," "the") that do not add value.

- Stemming: Reduces words to their root form (e.g., "coughing" becomes "cough").

Step 3: The normalized symptoms are matched with a predefined symptoms dictionary.

- Handling Spelling Mistakes: The system employs fuzzy matching to account for minor spelling errors, ensuring that similar-sounding or misspelled symptoms are correctly identified.

Step 4: The model predicts the disease based on the input vector created from the matched symptoms.

Step 5: receive a potential diagnosis and tailored recommendations

LAB REPORT ANALYSIS

The OCR Lab Report Analysis feature allows users to upload an image of their lab reports. The system processes the image to extract relevant text, summarize it, and highlight key information using advanced NLP techniques.

OCR Lab Report Analysis

Upload an image file to extract text using OCR.

Choose an image file...



Drag and drop file here

Limit 200MB per file • PNG, JPG, JPEG, BMP

Browse files



Home

Symptoms

Lab Report

X-Ray Report

Hiren Shah (PERSON)

Normachrami (ORG)

Norma (PERSON)

Mumber (PERSON)

Summary of Extracted Text:

Text Summary

"Peripheral BLOOD SMEAR EXAMINATION" was the subject of a blood test. The blood sample was collected at Siuvar Bungalow, 86 Road, 86 road, in Mumber. The sample was taken from a blood sample taken from the blood sample.

Read Summary


0:01 / 0:18



- **Image Upload:** Users can upload images in various formats (PNG, JPG, etc.) of their lab reports.
- **Image Preprocessing:** The uploaded image is converted to grayscale and binarized to enhance OCR accuracy.
- **Text Extraction:** Optical Character Recognition (OCR) is performed on the processed image to extract text content.(PYTESSERACT)
- **Named Entity Recognition (NER):** The extracted text is analyzed to identify and extract key information, such as dates, organizations, and quantities.(SPACY)
- **Text Summarization:** A summarization model, "[sshleifer/distilbart-cnn-12-6](#)," condenses the extracted text into a concise summary. This model is a distilled version of the BART (Bidirectional and Auto-Regressive Transformers) architecture, specifically fine-tuned for the task of summarization.
- **Text-to-Speech Option:** Users can listen to the summary via a text-to-speech feature, enhancing accessibility and user experience.(GTTS)

X-RAY ANALYSIS

The system classifies chest X-ray images to detect the presence of pneumonia, utilizing a trained Convolutional Neural Network (CNN) model. The user uploads an X-ray image, and the model returns a prediction along with a confidence score.



Step 1: User Uploads Chest X-ray Image

The user is prompted to upload a chest X-ray image in JPEG or PNG format.

Step 2: Model and Labels Loading

The pre-trained pneumonia classification model (xray_classifier.h5) is loaded, and class labels (e.g., "Normal", "Pneumonia") are read from a file.

Step 3: Image Preprocessing

The uploaded image is resized to (224, 224) pixels and normalized to fit the model's input requirements. This step ensures the image is in the right format for classification.

Step 4: Model Prediction

The preprocessed image is fed into the CNN model for classification. The model predicts whether the X-ray indicates "Pneumonia" or "Normal" along with a confidence score for the prediction.

Step 5: Displaying Results

The system displays the X-ray image along with the predicted class (Pneumonia or Normal) and the confidence score as a percentage.

Future Enhancements for BAYMAX

Improve Accuracy of CNN model ,Better User Interface and Optimizing the System


User-Specific Monitoring:

- Tailor health insights based on individual user data and preferences.
- Implement adaptive algorithms that adjust recommendations in real-time.


Integration of Llama 3:

- Utilize Llama 3 for enhanced natural language processing capabilities.
- Improve user interactions and data interpretation through advanced AI models.

 BAYMAX

 Home

 **Symptoms**

 Lab Report

 X-Ray Report

Symptom-based Medical Prediction

Enter Symptoms (comma-separated):

Predict

The background is a solid dark red color. In the four corners, there are decorative elements consisting of thin, light red lines that resemble circuit traces or a stylized tree structure, with small circles at the end of the lines.

THANK YOU

**PROJECT BY
RAJA PRASANNA M
PGA 22**